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Email info@medicooler.nl. You will get a personal answer from a specialist.

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The LifeinaBox is a portable mini fridge that actively keeps your medication between 2 and 8 degrees. At home, in your hotel, at the office or in the car: on mains power, via your car's 12V socket or on its built-in battery. Read this manual once from start to finish before first use, and after that you'll have it ready in minutes.

The LifeinaBox at a glance

Feature	LifeinaBox
Temperature	2–8°C (± 1°C)
Interior	180 x 65 x 40 mm (468 cm ³)
Power	110–240V AC (mains), 12V DC (car) or 5V USB-C
Built-in battery	approx. 6 hours
Max. ambient temperature	32°C

If your medication fits inside, the LifeinaBox keeps it at the right temperature. Suitable for insulin, adrenaline, growth hormone, arthritis medication and more, in prefilled pens, syringes and vials.

In the box

- LifeinaBox mini fridge
- Mains adapter 110–240V with Type A, Type I, Type C (EU), Type G plugs
- 5V USB-C cable
- Car cable for the 12V socket (cigarette lighter)

Ready to use in four steps

1 Fully charge the battery

Connect the LifeinaBox to mains power and set the O/I switch on the bottom to ON. Press and hold the button on the side for about 5 seconds until the logo on the front switches off. When the cooler is off, the red light next to the connector indicates the battery status: if it is blinking, it is charging; if it is full, it is solid red. Let it charge for at least 2 hours, or simply overnight. Then switch it on with the side button and check in the Lifeina app that the battery is full, or the red light indicator.

2 Position it and switch it on

Place the LifeinaBox lying flat on a level, firm surface, keeping the ventilation openings on the bottom and side clear. Switch it on with a short press of the side button, or unplug it and plug it back in.

3 Wait until it reaches temperature

After about 20 minutes the interior is stable between 2 and 8 degrees. You'll see it on the logo at the front: it starts to "breathe", slowly fading in and out.

4 Add your medication

Take your medication straight from the fridge and place it inside. Travelling with glass vials? Put some soft padding between them so they don't knock together and break.

Getting the best from your LifeinaBox

The LifeinaBox cools just like your fridge at home: it draws heat out of the interior and a fan releases it outside. That heat needs somewhere to go.

- Always use it lying flat on a level surface. Standing upright or on a soft, uneven surface, the ventilation is blocked and the device overheats.
- Always keep the ventilation openings on the bottom and side clear.
- Never use it switched on in an enclosed space without air circulation, such as a closed bag, suitcase, cupboard or desk drawer.
- Keep it out of direct sunlight and away from heat sources.
- Don't place it on the floor: the fan then sucks up dust, which will eventually clog the ventilation.
- Only open the LifeinaBox when you need your medication and close it again straight away.
- Do you have the stand? Use it for better airflow. The cooler sucks up less dust, which improves its lifespan, and it works quieter and more efficiently on the stand.

Good to know: the thermostat switches cooling on as soon as the interior reaches 8 degrees and off again at 2 degrees, so your medication never drops below 2 degrees. The LifeinaBox works up to an ambient temperature of 32 degrees. The warmer it is, the harder it has to work and the shorter the battery lasts. Please note: at ambient temperatures below 0 degrees, the LifeinaBox does not protect your medication from freezing.

Power and battery

You can power the LifeinaBox in three ways:

- Mains: with the supplied adapter (110–240V), usable worldwide.
- Car: with the supplied cable in the 12V socket (or the less powerful 5V USB-C cable).
- Battery: the built-in battery keeps it cool for about 6 hours, for example during a power cut or a layover. To use it, set the O/I switch on the bottom to ON.

Extra batteries (LifeinaPower) each give you about 6 more hours and can be charged with an ordinary phone cable.

Good to know: using the LifeinaBox at home on mains power? Keep the O/I switch on ON. This way, the battery will automatically take over if the power fails, and the battery will charge while the cooler is plugged in. Only set the switch to OFF if you are leaving the cooler plugged in for more than 5 days without use, as the battery will be disconnected and will not charge. Always fully charge it before you travel. Only use the supplied cables. A different or modified cable can damage the device, put the effectiveness of your medication at risk and voids the warranty.

Battery alarm

When the battery is down to 20 percent, you'll hear a beep. Briefly press the side button to silence the alarm, and connect the LifeinaBox to mains power or your car as soon as possible so cooling isn't interrupted.

Lights and signals

Signal	Meaning
Logo slowly "breathing"	The interior is between 2 and 8 degrees
Side light steady blue	Connected to your phone or computer via Bluetooth
Side light flashing red/blue	No Bluetooth connection
Beep	Battery down to 20 percent
2 short beeps	Fan error: check for obstruction or dust

Prefer no light, for example at night in your hotel room? Switch the battery on (O/I switch to ON) and quickly press the side button twice. The logo switches off.

The Lifeina app

With the Lifeina app (iOS and Android) you can check, among other things, whether the battery is fully charged.

- 1 Switch on Bluetooth on your phone.
- 2 Switch the cooler on with the switch to ON (the blue light flashes next to the port).
- 3 Go to the app.
- 4 Press "Find my LifeinaBox".
- 5 The connection is now established.

On the road by car

- Connect the LifeinaBox to the 12V socket with the supplied car cable.
- In the car, always use the LifeinaTravel transport bag and hang the LifeinaBox on the back of a seat. This keeps the ventilation openings clear.
- Never put it on the car floor: the fan sucks up dust from the floor and will eventually clog.
- Don't leave it in a parked car: the temperature there quickly rises far above 32 degrees.
- The LifeinaBox uses a maximum of 25 watts and won't drain a fully charged battery while the engine is running. Every car and battery is different: if you're unsure about the condition of your battery, have it checked by your garage.

Travelling, flying and security checks

- The LifeinaBox is designed for use at home, in hotels, at the office and in the car. For air travel, a passive cooling bag is far more practical: no power needed, lightweight and it slips easily into your hand luggage. Our cooling bags keep your medication between 2 and 8 degrees from 12 to 36 hours depending on the product. You'll find all our cooling bags at www.medicooler.nl.
- Always keep your medication with you in the cabin, never in the hold. If your luggage goes missing or your flight is delayed, you still have your medication with you.
- Taking the LifeinaBox along as a fridge for your destination? Switch it off completely before packing it. Switched on in a closed bag or suitcase, it overheats.

- Want to use it during the flight? That's possible with the LifeinaTravel bag, which keeps the ventilation openings clear. Please note: most airlines do not allow the use of power banks, so you are limited to the maximum 6-hour battery life. You can charge it using the power outlet under your seat or in the headrest in front of you. Do not place the cooler in the overhead bin, as this will lead to suffocation and overheating.
- Under aviation rules, spare batteries such as powerbanks may only be carried in hand luggage, never in the hold. Check your airline's rules before you fly.
- Your medication is medically necessary and falls under the exemption to the 100ml rule. Show it separately at security if asked.
- Customs sometimes ask for a doctor's certificate for the medication itself. Travelling outside Europe? Check the destination country's rules on bringing medication in advance.

At your destination

Place the LifeinaBox lying flat on a level surface with clear ventilation, plug it into mains power and wait about 20 minutes until it reaches temperature. Then add your medication. This gives you your own medicine fridge in every hotel, so you don't depend on a minibar or a shared fridge. Need the battery the next day? Leave the O/I switch on ON so it fully charges overnight.

Care and storage

- Clean the LifeinaBox with a cloth lightly moistened with water or a mild detergent.
- Never submerge it in water and don't expose it to excessive moisture. The LifeinaBox is not waterproof.
- Don't drop it, and never open or modify it yourself.
- Don't bend or twist the power cables.
- Keep the ventilation openings clean and free of dust.
- Store it in a dry place, away from sunlight and dust. Fully charge the battery before your next trip.

Technical specifications

Specification	LifeinaBox
Outer dimensions	220 x 107 x 110 mm (W x D x H)
Interior	180 x 65 x 40 mm (468 cm ³)
Weight	1.4 kg
Voltage	12V DC (car) or 110–240V AC
Power consumption	max. 25 watts
Temperature range	2–8°C, ± 1°C
Built-in battery	approx. 6 hours
Insulation	CFC-free polyurethane foam
Cooling system	Wear-free thermoelectric cooling with fan
Material	ABS

Specifications are subject to change without notice due to technical improvements.

| Good to know: what affects the temperature

The LifeinaBox has been designed and tested with the greatest care. Even so, some conditions affect the performance of any portable fridge and are beyond our control. Please keep these in mind:

- the ambient temperature (max. 32°C)
- clear ventilation on the bottom and side
- a level, horizontal position
- exposure to direct sunlight
- how often you open the LifeinaBox
- the temperature of your medication when you put it in
- the battery charge level
- dust in the ventilation openings

DISCLAIMER

Responsibility for storing medication at the correct temperature lies entirely with the user. MediCooler shall under no circumstances be liable for any direct or indirect damage to property, medication or health resulting from the use or misuse of the LifeinaBox, its batteries or cables. Having doubts about your purchase after receiving it? Please contact our customer service at info@medicooler.nl before first use to return the LifeinaBox.

| Service and warranty

In addition to your statutory rights, Lifeina provides a one-year manufacturer's warranty from the date of purchase covering internal faults and manufacturing defects. The warranty does not cover damage caused by dropping or impact, incorrect handling, third-party intervention, natural disasters or water. Questions about your trip, your medication or your LifeinaBox? Email info@medicooler.nl. You'll get a personal reply from a specialist.

We wish you a good trip. Your medication stays at the right temperature.

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